

**INSTRUCTIONS**  
for  
**Electric Penis Pump**  
**A Male Enhancement Vacuum System**  
For  
**Penis Enlargement**

The results are faster than using a hand pump as the electric pump creates a more uniform and a uniform and constant pressure not available with any hand pump. This maximizes penis growth.

Your Electric Penis Pump uses a state-of-the-art vacuum system to create the most reliable and efficient penis enhancement system.

- A Two-Inch Diameter Expansion Cylinder
- An Electric Pump for Constant Pressure
- An Air Pressure Gauge
- A Safety Release Valve
- Three Rubber Constriction Rings
- A Leather Constriction Ring for Post Therapy Erections

**INTRODUCTION**

Thank you for purchasing the Electric Penis Pump. You are now on your way to regaining self-confidence and enhancing your intimate life. This system was specifically designed for assisting in Penis Enlargement.

**CRITICAL**

Please do not to use this pump until you slowly read all the instructions. Thank you.

**BEFORE PUMPING**

Familiarize yourself with the Electric Penis Pump, and its' operation.

Check out the ease of using the safety release valve.

Do Not Ejaculate for several hours prior to or during a session or even several hours afterwards. This allows maximum expansion without swelling the penis.


Massage and gently stretch the penis and scrotum for 3 to 5 minutes.

Generously lubricate the head. Shaft and base of the penis. For a secure seal of the cylinder to the body, apply a lubricant directly on the pubic hair.

Having a visual stimulant handy or a creative imagination is helpful. A friend to assist you with visual stimulation or even light touching is also very helpful.

## OPERATION OF THE ELECTRIC PENIS PUMP

Familiarize yourself with the pump, and how to use it. More importantly, you need to know the location and how to use the quick-release valve.



Connect the cylinder to the Electric pump by simply inserting the long, clear tubing onto the nozzle sticking out of the pump. This not easy to do unless you follow these directions. Always soften the tubing before you put it onto the nozzle. See picture above. The tubing is easily softened by putting it into very hot or even boiling water for a minute or two. Shake off any excess water, and then you can easily insert the tubing. Any remaining water will dry rapidly, and the tube will shrink tightly to provide a leak-free seal.

Slide one constriction ring over the cylinder, and then down to about 1/4 inch ( about 7 mm) from the base of the cylinder. You have received 3 constriction rings. Try the largest ring first---use the one that is the most comfortable.

Prior to beginning pumping, stimulate the penis with your hand to achieve your best possible erection. Use a light amount of lubricant as it will assist stimulation. It will also inhibit the penis from sticking in the cylinder as it is being pumped.

Place the cylinder over the penis, then slide the constriction ring off of the cylinder and onto the penis.

### **CRITICAL----Please Read Slowly**

**The Electric Penis Pump is almost completely assembled when you receive it. The release valve handle needs to be turned totally to the open position, which means you gently turn it counter-clockwise until it reaches the end point and stops turning. Handle it gently; it will turn very easily. Do not turn it roughly or past the the end point. Now insert your penis into the cylinder, and slowly turn the release value clockwise with the pump operating. Turn it clockwise very slowly, until you are able to change your flaccid penis into an erect penis.**

**After you have completed pumping, and you have a fully erect penis, slide the ring off the cylinder and onto the base of the shaft of the penis while simultaneously releasing the pressure from the cylinder and removing the cylinder from the penis. The ring will then help you to stay erect and engage in sex.**

**Subsequently, wash penis and groin area thoroughly with an antibacterial soap and warm water.**

**DO NOT keep constriction ring on for more then 45 minutes.**

**DO NOT sleep with constriction ring on. If penis becomes bruised as a result of using the constriction rings, discontinue usage until bruising heals. If there is any pain, consult a physician.**

**You can further increase the enlargement process rate by the use of the leather strap. Secure snugly but comfortably the leather penis strap under the scrotum and snapped over the penis after your pumping session. Erection constriction rings should never be worn whenever the strap is being used.**

**The gauge is a technical tool that measures pressure, and the dial should never be at any value that is uncomfortable for you to use. If you look at it carefully, you will see large numbers printed on the dial face, and many small increments in between the numbers. As you begin your initial vacuum pumping, always use a pressure that is very comfortable. In fact, it should always be comfortable. During the first month or so, always use the same pressure. After about one month, you can increase the pressure. Only increase the level very slightly. For example, where ever the needle of the gauge is pointing during the first month, it only needs to move one to three small increments on the dial. Do not try to increase from one digit to another. Increasing pressure must be done very gradually.**

**You should be pumping for about 5 minutes per session in the first 6 week period,**

**and subsequently after this time, you can do a second 5 minutes pumping session. It is best to use the Electric Pump three days within every 7 to 8 days, skipping one or two days in between. Never use the penis cylinder two days in a row--- Like body building , it is important to allow the tissues to heal.**

**You should see some results after about 6 weeks or so, and permanent results after 3 to 4 months,**

### **CAUTION**

**If you have symptoms of erectile dysfunction (i.e., inability to achieve an erection sufficient for sexual intercourse), consult your physician prior to using this device to avoid a potentially harmful delay in diagnosing any of the most common causes of this condition, such as diabetes, multiple sclerosis, cirrhosis of the liver, chronic renal failure, or alcoholism.**

**Apply the minimum amount of pressure necessary to achieve an erection. Stop using if any pain occurs. Do not use if under the influence of alcohol or drugs. Misuse may result in aggravation of existing medical conditions, swelling, of the penis, or permanent injury to the penis. Do not use electrically powered pumps in or near water. Consult a doctor and discontinue use if you experience prolonged pain, redness, irritation or bruising.**

#### **Disclaimer:**

**This item is sold as a novelty item only. No medical claims are warranted or implied by the use of this product.**